

## cDNA library production by OTTR using low-cycle PCR for indexing

*Use pipettes designed for low-volume measurements (e.g. a p2). If low-volume pipettes are not available, please double all volumes listed below for greater accuracy.*

*Sequential buffer additions to the OTTR reaction tube may be done at room temperature or on ice; ice should be used when performing many reactions in parallel.*

### STEP 1 Label template 3' ends using ddRTPs.

- a. Adjust input single-stranded nucleic acid to **9  $\mu$ L** using sterile water.\*
- b. Add **2  $\mu$ L** Buffer 1A.
- c. Add **1  $\mu$ L** Buffer 1B.
- d. Add **1  $\mu$ L** Tail Pol.
- e. Lightly mix and quick-spin.
- f. Incubate at **30°C** 90 min.
- g. Add **1  $\mu$ L** Buffer 1C.
- h. Lightly mix and quickly spin.
- i. Incubate at **30°C** 30 min.
- j. Incubate at **65°C** 3-5 min, then place on ice.

*Safe to stop; store at -20°C.*

### STEP 2 Eliminate unincorporated ddRTPs.

- k. Add **1  $\mu$ L** Buffer 2
- l. Add **1  $\mu$ L** diluted rSAP
- m. Incubate at **37°C** 15 min, then place on ice.

### STEP 3 Heat-neutralize phosphatase.

- n. Add **1  $\mu$ L** Buffer 3.
- o. Lightly mix and quick-spin.
- p. Incubate at **65°C** 5 min, then place on ice.

*Safe to stop; store at -20 °C.*

### STEP 4 cDNA synthesis.

- q. Add **1  $\mu$ L** Buffer 4A.\*\*
- r. Add **1  $\mu$ L** Buffer 4B.
- s. Briefly vortex then quick-spin to ensure good mixing of viscous buffers.
- t. Add **1  $\mu$ L** Relay Pol.
- u. Lightly mix and quick-spin. *Place on ice if not rapidly transferred to 37°C.*
- v. Incubate at **37°C** 15 min.
- w. Incubate at **65°C** 5 min, then cool.

\* *Input RNA should not be strongly buffered and should lack >0.1 mM chelating agent (e.g. EDTA)*

\*\* *Buffer 4A is viscous; check that the entire tube contents is collected & well mixed at the bottom of the kit buffer tube, and watch your pipet tip to not bring extra volume into OTTR reactions.*

**EXAMPLE OF OPTIONAL cDNA CLEAN UP      Reagents not provided**

- a. Add 0.5  $\mu$ L RNase H (NEB Thermostable) + 0.5  $\mu$ L RNase A (Sigma).
- b. Incubate at 50°C ~10 min.
- c. Add 80  $\mu$ L stop mix (50 mM Tris pH 8, 20 mM EDTA, 0.2% SDS).
- d. Extract with equal volume phenol:chloroform:isoamyl alcohol (25:24:1).
- e. Move the upper aqueous phase to a new tube containing 10  $\mu$ L 3 M NaOAc pH 5.2 and 20  $\mu$ g glycogen. Mix. Add 350  $\mu$ L 100% ethanol. Mix.
- f. Chill in dry ice bath 2-3 min. Spin 15 min at max microfuge speed. Remove supernatant and wash pellet with 1 mL 70% ethanol. Spin again.
- g. Remove all supernatant and dry pellet for 5 min. Resuspend for next use.

**EXAMPLE OF PCR INDEXING FOR ILLUMINA NGS      Reagents not provided**

For 50  $\mu$ L PCR reactions, make a master mix combining for each reaction:

- 16.5  $\mu$ L sterile water
- 10  $\mu$ L 5X Q5 buffer (NEB)
- 1  $\mu$ L 10 mM dNTP mix
- 1  $\mu$ L 10  $\mu$ M P5(i5 index)R1 primer
- 1  $\mu$ L 10  $\mu$ M P7(i7 index)R2 primer (complement strand)
- 0.5  $\mu$ L Q5 polymerase (NEB)

Add 30  $\mu$ L master mix above to 20  $\mu$ L cDNA + sterile water

Amplify in PCR machine with heated lid:

98°C, 1 min

8X Cycle (increase cycles for very low input)

98°C, 20 sec

65°C, 20 sec

72°C, 10 sec

72°C, 5 min



**ILLUMINA NGS READ ANALYSIS**

*Read 1 orientation is input 5' to 3' with the first 7 nucleotides giving the adaptor template unique molecular identifier (UMI).*

*Read 2 orientation is input 3' to 5' with the first nucleotide from the adaptor primer 3' end.*

## BUFFER COMPONENTS

Buffer 1A	140 mM Tris-HCl (pH 7.5), 1 M KCl, 14 mM DTT, 35% PEG-8000
Buffer 1B	28 mM MnCl <sub>2</sub> , 3.5 mM ddATP
Buffer 1C	3.5 mM ddGTP
Buffer 2	80 mM MgCl <sub>2</sub>
Buffer 3	100 mM EGTA
Buffer 4A	10 mM MgCl <sub>2</sub> , 900 mM KCl, 40% PEG-6000; 4 mM dGTP, 0.8 mM dTTP, 0.8 mM dCTP, 0.04 mM dATP, 3 mM DAP-TP
Buffer 4B	Adaptor oligonucleotide mix: 1.8 μM annealed primer duplex mix <b>U5PD</b> , 3.6 μM adaptor template <b>U3T</b>

## ADAPTOR OLIGONUCLEOTIDES

U5PD oligos	<u>c5pUNI_T</u> (primer) 5Cy5/GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
	<u>c5pUNI_C</u> (primer) 5Cy5/GTGACTGGAGTTCAGACGTGTGCTCTTCCGATCC
	<u>c5pc</u> (primer complement strand) rGrArUrCrGrGrArArGrAmGmCmAmCmAmCmGmUmCmUmG mAmAmCmUmCmCmAmGmU/3SpC3/
U3T oligo	<u>c3tUNI-UMI*</u> (cDNA 3' adaptor template) 5Cy5/ACACTCTTCCCTACACGACGCTCTTCCGATCTrNrNrNr NrNrYrC

## ENZYMES

Tail Pol	BoMoC(ed) F753A 10 μM AddGene Plasmid #185710
Relay Pol	BoMoC(ed) 10 μM AddGene Plasmid #186461
rSAP	NEB M0371S diluted 1 vol: 1 vol into 50% glycerol, final 0.5 U/μL

## For reference:

Upton HE, Ferguson L, Temoche-Diaz MM, Liu XM, Pimentel SC, Ingolia NT, Schekman R, Collins K. *Low-bias ncRNA libraries using ordered two-template relay: Serial template jumping by a modified retroelement reverse transcriptase. PNAS* 118: e2107900118 (2021).

Pimentel SC, Upton HE, Collins K. *Separable structural requirements for cDNA synthesis, nontemplated extension, and template jumping by a non-LTR retroelement reverse transcriptase. J Biol Chem* 298: 101624 (2022).